

Please note that the following setup is calculated for 200 mL inoculated medium that gives 50 mL of Taq solution.

Day 1

- inoculate 3 ml LB containing Carbenicillin (25 µg/mL) with glycerol stock
- grow overnight or from 8 am till 3 pm at 37°C

Day 2

- take 2 ml of the preculture, inoculate 200 ml LB/Carb in a 2 L Erlenmeyer flask or scale up appropriately
- immediately prepare a fresh glycerol stock culture from the remaining mL for further use
- let grow till culture reaches OD₆₀₀ to 0.3 (approx. 3 hrs)
- add 0.5 mM IPTG (1 ml/200 mL culture 100 mM IPTG)
- let the culture grow o/n @ 37°C (16 h)

Day 3

- harvest the cells by centrifugation 10 min @ 4000 rpm and discard the supernatant
- store the pellet at -20°C (essential step!), at least 1 h or overnight to break the cells

Day 4

- resuspend the cells in buffer A: 6 ml/200 mL culture
(OVERALL VOL. is now approx. 6 mL/200 mL culture)
- sonify suspension on ice for 2 min @ cycle 5 and 20% efficiency
- add 10 mM MgCl₂ (60 µl/200 mL culture 1 M MgCl₂)
- add DNase I. For 200 mL culture, we used one tenth of a complete vial RNase-Free DNase Set (Qiagen, 1500 Units, Cat. No. 79254, 65,00 €, take up in 550 µL of the provided RNase-free water) or 0.2 µl/ml DNase (unknown provider/concentration)
- incubate 30 min at 37°C in shaking sterilized waterbath or on the roller (suspension gets yellowish)
- add buffer B: 6 ml/200 mL culture
(OVERALL VOL. is now approx. 12 mL/200 mL culture)
- heat to 75°C for 60 min in a shaking sterilized water bath until precipitation occurs
- centrifuge at 13.000 rpm @ 4°C for 20 min
- mix the yellowish supernatant with one volume of storage buffer 50% glycerol: approx. 12 mL/200 mL culture
(OVERALL VOL. is now approx. 24 mL/200 mL culture)
- mix now with one volume of the cloudy storage buffer 75% glycerol approx. 24 mL/200 mL culture (2 times the initial supernatant volume)
(OVERALL VOL. is now approx. 48 mL/200 mL culture)
- aliquot to 500 µL and store the Taq solution (final glycerol conc. is 50%) at -20°C

Day 5

- check the quality by using different concentrations of the freshly prepared Taq compared to the previous stock

Normally 0.3-0.5µl of this preparation in a 50µl PCR reaction is sufficient. High amounts of *Taq* are inhibitory

Buffer A

10 ml (needed for 200 mL culture are 6 mL)

		<u>check here:</u>
50 mM Tris/HCL pH 7.9	1 M Tris pH 7.9 → 500 µl	<input type="checkbox"/>
50 mM Glucose (M = 180,16 g/mol)	→ 90,8 mg	<input type="checkbox"/>
1 mM EDTA	0,5 M EDTA → 20 µl	<input type="checkbox"/>
4 mg/ml Lysozyme	→ 40 mg	<input type="checkbox"/>
	ad 10 ml H ₂ O	<input type="checkbox"/>

Buffer B

10 ml (needed for 200 mL culture are 6 mL)

		<u>check here:</u>
10 mM Tris/HCL pH 7.9	1 M Tris pH 7.9 → 100 µl	<input type="checkbox"/>
50 mM KCL	2,5 M KCL → 200 µl	<input type="checkbox"/>
1 mM EDTA	0,5 M EDTA → 20 µl	<input type="checkbox"/>
0,5% Nonidet P-40	10 % Nonidet → 500 µl	<input type="checkbox"/>
0,5% Tween 20	10 % Tween 20 → 500 µl	<input type="checkbox"/>
	ad 10 ml H ₂ O	<input type="checkbox"/>

10X Storage buffer

(used to prepare the two following glycerol containing buffers)

10 ml (needed for 200 mL culture are 5 mL)

final conc. in 1X storage buffer (not in the 10X stock!):

		<u>check here:</u>
50 mM Tris/HCL pH 8	1 M Tris pH 8 → 2.5 ml	<input type="checkbox"/>
100 mM NaCl	5 M NaCl → 1 ml	<input type="checkbox"/>
1 mM EDTA	0,5 M EDTA → 100 µl	<input type="checkbox"/>
0,5 mM DTT (M:154,25 g/mol)	50 mM DTT → 500 µl	<input type="checkbox"/>
1% Triton X-100	100% Triton X-10 → 500 µl	<input type="checkbox"/>
	ad 5 ml H ₂ O	<input type="checkbox"/>

50% Glycerol storage buffer

20 mL (needed for 200 mL culture are 12 mL)

	<u>check here:</u>
12,5 g Glycerol (= 10 ml; density = 1,249 g/mL)	<input type="checkbox"/>
2 ml 10X storage buffer	<input type="checkbox"/>
ad 20 ml H ₂ O	<input type="checkbox"/>

75% Glycerol storage buffer (*sometimes very cloudy appearance*)

30 mL (needed for 200 mL culture are 24 mL)

	<u>check here:</u>
28,125 g Glycerol (= 22,5 ml; density = 1,249 g/mL)	<input type="checkbox"/>
3 ml storage buffer (10x)	<input type="checkbox"/>
ad 30 ml H ₂ O	<input type="checkbox"/>